

LIFEPAK® 12 DEFIBRILLATOR/MONITOR

Works like
you work.™



PHYSIO
CONTROL



u make a difference.

Your job demands more.

Use a tool that can tackle today's patient care needs and adapt to tomorrow's challenges.

The gold standard for more than 30 years, LIFEPAK 12 defibrillators are continually evolving to keep pace with the patient care. The LIFEPAK 12 defibrillator/monitor is a parameter therapeutic and diagnostic function portable device.

Over 80,000 LIFEPAK 12 defibrillator/monitors are on rescue rigs and in hospitals worldwide. Feedback from the global community keeps us innovating — adding new features you in your lifesaving work.



Trusted

EVOLUTION

Keeping pace as patient care evolves.



Here are highlights of advances since its first release:

Defibrillator/monitor
for acute cardiac care
provided diagnostic
capabilities.
It also does the 12.

1998

The LIFEPAK 12 defibrillator/monitor revolutionizes acute cardiac care, with expanded diagnostic and monitoring capabilities.

1999

LIFEPAK 12 defibrillator/monitor is enhanced with ADAPTIV™ biphasic technology up to 360J, NIBP and CO₂ monitoring capabilities.

2001

Among many features added to the 12 are ST Monitoring up to 8 hours and invasive pressure monitoring.

2004

Bluetooth® capabilities enable wireless transmission of 12-lead ECGs.

2007

cprMAX™ technology provides increased flexibility for protocols to maximize CPR. STEMI Management technology enables secure and flexible flow of ECG data, linking EMS and hospitals for improved STEMI Treatment.



Complemented by a rich range of services and options

Technical Field Service

We offer one of the largest and best-trained networks of technicians in the industry. On call 24 hours a day, 7 days a week (North America), your phone call within two hours, to work with you to quickly assess the best solution. An integral part of your team, our field service technicians have a tenure of more than 12 years and log an average of 2,400 field hours with you to design a customized service offering that meets your needs. We offer a comprehensive plan covering repair, preventative maintenance and limited to repair-only or inspection-only.

Training

Whether you are taking delivery of your first LIFEPAK 12 defibrillator or new options, your sales rep will provide inservice training to help you understand your Physio-Control products. Specialized training — ranging from live webcasts to on-site classes — is also available for features such as monitoring and 12-lead ECG. Continuing education credits are available.

Grant Consultation

This complimentary program helps Physio-Control customers streamline the grant and writing process. Free assistance is provided to communities, fire and EMS agencies in applying for grants to upgrade equipment, purchase programs, and provide training. Customers receive guidance on identifying funding sources, timing grant applications, and crafting a grant proposal to improve

Accessories

demands of your changing job.

use when cond counts

e you work—in the most demanding

is you quickly scan clutter-free waveforms

drop test/impact test of 18 inches;

N MONITORING

lead ECG acquisition and transmission, advance the state of the art, with ST changes in a patient's 12-lead ECG.

cor/monitor on the market today* with because ECGs (and the diagnosis) can quickly, the device takes a series of and alerts you to changes in a patient's

status breath by breath with patented technology and FilterLine® accessories in high humidity. EtCO₂ monitoring is and nonintubated patients.

is allows for evaluation of changes in nt response to therapy over time.

ity offers accurate and stable oxygen



MANUAL DEFIBRILLATION

1. Push ON. Apply conductive gel to hand paddles or push ANALYZE. Push CHARGE. Push SHOCK to deliver energy.


AED OPERATION

- Push ON.
- Push ANALYZE.
- Push SHOCK when directed to deliver energy.

PACER OPERATION

- Push PACER to turn pacer on.
- Push RATE to set rate and adjust up or down as needed.
- Push CURRENT to set current and adjust up or down as needed.

Medtronic
LIFEPAK 12
BIPHASIC DEFIBRILLATOR/MONITOR



1 ON

2 ADVISORY ENERGY SELECT

3 ANALYZE CHARGE SHOCK

LEAD SIZE SYNC

NIBP PACER

ALARMS RATE

OPTIONS CURRENT

EVENT PAUSE

Home Screen

12-LEAD

TRANSMIT

CODE SUMMARY

PRINT

CO2

SP02

NIBP

ECG

P1

P2



standard. Raising the bar.

Improving treatment with secure flow of ECG data

Mortality for Acute Coronary Syndrome (ACS) has been shown to increase 40% if door-to-balloon time stretches from 90 minutes to 120 minutes.⁴

Transmitting 12-lead ECGs from the field so hospitals can confirm STEMI diagnosis and activate the cardiac cath lab before the patient arrives can help you meet the ACC/AHA 90-minute door-to-balloon guideline for STEMI patients.⁵ The LIFENET STEMI Management Solution links EMS response with hospitals via the internet to improve treatment and management of patients with ACS. Using cutting edge technology, 12-lead ECG data is sent from the 12 to a network of receiving targets in the ED, cath lab, or anywhere the LIFENET Client software is installed, so patient data can be shared quickly. This technology enables hospitals to make patient care decisions such as referral to a PCI capable hospital, or cath lab activation, all while you are still en route.

Studies have reported a significant association between prehospital 12-lead ECGs and shorter door-to-balloon times. Two recent studies found the effect was strongest when the cath lab was activated while the patient was still en route to the hospital.^{6,7}



LIFENET® STEMI Management Solution technology facilitates a seamless, secure flow of ECG data between prehospital, emergency, and hospital treatment centers, enabling you to improve patient outcomes and reduce false-positive cath lab activations. The LIFENET STEMI care network linking EMS responders with your gateway devices like SmartPhone or tablet, in the field focus on patient care, a gateway to diagnostic quality ECGs wirelessly from the field to the proper destination, helping you meet the 90-minute door-to-balloon guideline for STEMI. The system works on a variety of wireless devices and requires no dedicated hardware. Applications are a state-of-the-art datacenter that maximize data availability—are provided on a subscription basis—support are available to meet your needs.

CODE-STAT™ 7.0 Data Review Software **Advanced CPR Analytics.** This powerful software annotates chest compressions onto the ECG report and calculates CPR statistics according to 2005 AHA Guidelines. The software simplifies reporting by consolidating all display data and outcome data into a single e-file. With download, review, manage and analyze data from multiple LIFEPAK defibrillators, facilitates quality analysis and business creation of benchmarking and trending system's performance.

DT EXPRESS™ Data Transfer Software Windows-based software application for LIFEPAK devices. The software makes critical event and waveform data to your patient data, print out a hardcopy report on disk. For storage and on-screen view to CODE-STAT 7.0 data review software.

Powerful

Series has five main features available continuously	Operating capability for				NICd (FASTPAK® FASTPAK NiCd battery) or			
	tomize the device	ecute device	waveforms for					
with batteries (unit and lbs) Add 0.3 lbs	power (<30				with			
	ected							
Low battery icon at message in status area is indicated, device when both batteries are a voice prompt	power (<30				with			
	ected							

Standard Paddles (hard): 0.9kg (1.9 lbs)	
Height: 31.7cm (12.5 in)	
Width: 39.6cm (15.6 in)	
Depth: 23.1cm (9.1 in)	
DISPLAY	
Size (active viewing area):	
LCD: 140.8mm (5.5 in) wide x 105.6mm (4.2 in) high	
EL: 165.1mm (6.5 in) wide x 123.8mm (4.9 in) high	
Resolution:	
640 x 480 black and white LCD	
640 x 480 amber and black EL display	
User selectable LCD contrast	
Displays a minimum of 4 seconds of ECG and alphanumeric for values, device instructions or prompts	
Option to display one or two additional waveforms	
Waveform Display Sweep Speed: 25mm/sec for ECG and 12.5mm/sec of CO2	

DATA MANAGEMENT	
The device captures and stores patient data, events (including waveforms and annotations), user test results and continuous ECG waveform records in internal memory.	
The user can select and print reports and transfer the stored information via an internal modem through landline or mobile phones.	
Report Types: Three format types of CODE SUMMARY™ critical event record (short, medium and long)	
- Initial ECG (except short format) - Automatic capture of vital signs measurements every 5 minutes 3-channel or 4-channel 12-lead ECG report - Continuous waveform records (transfer only) - Trend Summary – Includes patient information, vital signs log and vital signs graphs - Vital Signs – Includes patient information, event and vital signs log - Snapshot – Includes patient information and 8 seconds of ECG captured at the time of transmission	
Memory Capacity: Two full-capacity patient records that include:	
CODE SUMMARY critical event record – up to 100 single waveform events	
Continuous Waveform – 45-minute continuous ECG record	

COMMUNICATIONS	
The device is capable of transferring data records by internal modem, external EIA/TIA modem, cellular modem or serial connection	
Bluetooth wireless data transfer to cell phone to LIFENET RS receiving station	
Supports EIA/TIA-602 compatible modems using Xon/Xoff or RTS/CTS flow control at 9600 to 38400 bps	
EIA/TIA-RS232E compatible at 9600, 19200, 38400 and 57600 bps	
Group III, Class 2 or 2.0 fax	

ECG	
ECG is monitored via several cable arrangements:	
A 3-wire cable is used for 3-lead ECG monitoring	
A 5-wire cable is used for 7-lead monitoring	
A 10-wire cable is used for 12-lead acquisition. When the chest electrodes are removed, the 10-wire cable functions as a 4-wire cable.	
Standard paddles or QUIK-COMBO pacing/defibrillation/ECG electrodes or FAST-PATCH® disposable defibrillation/ECG electrodes are used for paddles lead monitoring	
Lead Selection: Leads I, II, III, (3-wire ECG cable)	
Leads I, II, III, AVR, AVL and AVF acquired simultaneously (4-wire ECG cable)	
Leads I, II, III, AVR, AVL, AVF, V1 (Labeled “C” on 5-wire ECG cable)	
Leads I, II, III, AVR, AVL, V1, V2, V3, V4, V5 and V6 acquired simultaneously, (10-wire ECG cable)	
ECG Size: 4, 3, 2.5, 2, 1.5, 1, 0.5, 0.25 cm/mV (fixed at 1 cm/mV for 12-lead)	
Heart Rate Display: 20 to 300 bpm digital display	
Out of Range Indication: Display symbol “- - -”	
Heart symbol flashes for each QRS detection	
Continuous Patient Surveillance System (CPSS): In advisory mode while Shock Advisory System™ is not active, CPSS monitors the patient, via paddles or Lead II ECG, for potentially shockable rhythms	
Analog ECG Output: 1V/mV x 1.0 gain	
Common Mode Rejection: 90dB at 50/60Hz	

SpO2	
MASIMO SET Sensors	
Saturation Range: 1 to 100%	
Saturation Accuracy: (70–100%) (0–69% unspecified)	
Adults/Pediatrics:	
+/- 2 digits (during no motion conditions)	
+/- 3 digits (during motion conditions)	
Neonates:	
+/- 3 digits (during no motion conditions)	
+/- 3 digits (during motion conditions)	
Dynamic signal strength bar graph	
Pulse tone at the onset of the pleth waveform	
SpO2 Update Averaging Rate: User selectable 4, 8, 12 or 16 seconds	
SpO2 Measurement: Functional SpO2 values are displayed and stored	
Pulse Rate Range: 25 to 240 pulses per minute	
Pulse Rate Accuracy: (Adults/Pediatrics/Neonates)	
+/- 3 digits (during no motion conditions)	
+/- 5 digits (during motion conditions)	
SpO2 waveform with autogain control	

NIBP	
Oscillometric measurement	
Systolic Pressure Range: 30 to 245mmHg	
Diastolic Pressure Range: 12 to 210mmHg	
Units: mmHg, kPa	
Mean Arterial Pressure Range: 20 to 225mmHg	
Blood Pressure Accuracy: maximum mean error of ± 5mmHg with a standard deviation no greater than ± 8mmHg	
Pulse Rate Range: 30 to 200 pulses per minute	
Pulse Rate Accuracy: 2 pulses per minute or 2%	

EtCO2	
Microstream technology	
Measurement range: 0 to 99mmHg	
Display: CO2 waveform and EtCO2 numerics	
Units: mmHg, kPa, %; user selectable	
Automatic ambient pressure compensation	
CO2 Accuracy (>20 minutes): 0 to 38mmHg; ± 2mmHg@9 to 99mmHg; ± 5% of reading + 0.08% for every 1mmHg	
Warm Up Time: 30 seconds (typical), 180 seconds max	
Response Time: 2.9 seconds (includes delay time and rise time)	
Respiration Rate Range: 0 to 60 breaths per minute	
Respiration Rate Accuracy: 0 to 40 bpm: ± 1 bpm, 41 to 60 bpm: ± 2 bpm	

Invasive Pressure (2 channels)	
Measurement Range: -30 to +300mmHg in six user selectable ranges	
Display: IP waveform and numerics	
Units: mmHg, kPa	
User-selectable Labels: ART, PA, CVP, ICP, LAP	
Transducer Type: Strain-gauge resistive bridge	
Transducer Sensitivity: 5mV/V/mmHg	
Bandwidth: 0 - 30 Hz (<3dB)	
Numeric Accuracy: ± 1mmHg or 2% of reading, whichever is greater, plus transducer error	
Leakage Current: Meets ANSI/AAMI/IEC requirements	
Trend	
Display: Choice of HR, SpO2(%), EtCO2, RR, NIBP, P1, P2, ST shown in channels 2 or 3	
Time Scale: Auto, 30 minutes, 1, 2, 4 or 8 hours	
Duration: Up to 8 hours with -06 Memory PCB or later. Reduced storage capacity with earlier versions.	
ST Segment: After initial 12-lead ECG analysis, automatically selects and trends lead with the greatest ST displacement	

ALARMS	
Quick Set: Activates alarms for all parameters	
VF/VT Alarm: Activates continuous CPSS monitoring in Manual Mode	
Apnea Alarm: Occurs when 30 seconds have elapsed since last detected respiration	

INTERPRETIVE ALGORITHMS	
12-Lead Interpretive Algorithm: GE Medical 12SL, Includes AMI statement	
PRINTER	
Prints continuous strip of the displayed patient information	
Paper Size: 50mm (2.0 in) or optional 100mm (3.9 in)	
Print Speed: 25mm/Sec +/- 5% (measured in accordance with AAMI EC-11, 4.2.5.2)	
Delay: 8 seconds	
Autoprint: Waveform events print automatically (user configurable)	
Optional 50mm/sec timebase for 12-lead ECG reports	

FREQUENCY RESPONSE	
Diagnostic: 0.05 to 150Hz or 0.05 to 40Hz (user configurable)	
Monitor: 0.67 to 40Hz or 1 to 30Hz (user configurable)	
Paddles: 2.5 to 30Hz	
Analog ECG Output: 0.67 to 32Hz (except 2.5 to 30Hz for Paddles ECG and 1.3 to 23Hz for 1 to 30Hz monitor frequency response)	

DEFIBRILLATOR	
Waveform: Biphasic truncated exponential with voltage and duration compensation for patient impedance	
Energy Accuracy: ±1 joule or 10% of setting, whichever is greater, into 50 ohms	
±1 joule or ±5%, whichever is greater, of 50 ohm value into 25 to 200 ohms*	
Paddle Options: QUIK-COMBO pacing/defibrillation/ECG electrodes (standard)	
FAST-PATCH disposable defibrillation/ECG electrodes (optional)	
Standard Paddles (optional)	
Internal Handles with discharge control (optional)	
External Sterilizable Paddles (optional)	
Cable Length: 2.4m (8 ft) long QUIK-COMBO cable (not including electrode assembly)	

Manual	
Energy Select (Biphasic): 2, 3, 4, 5, 6, 7, 8, 9, 10, 15, 20, 30, 50, 70, 100, 125, 150, 175, 200, 225, 250, 275, 300, 325, and 360 joules or user configurable sequence 100 to 360 joules	
Energy Select (Internal): 2, 3, 4, 5, 6, 7, 8, 9, 10, 15, 20, 30 and 50 joules	
Charge Time: Charge time to 360J in less than 10 seconds, typical	
Synchronous Cardioversion: Energy transfer begins within 60ms of the QRS peak	

AED	
Shock Advisory System (SAS): an ECG analysis system that advises the operator if the algorithm detects a shockable or non-shockable ECG rhythm. SAS acquires ECG via therapy electrodes only.	
Shock Ready Time: Using a fully charged battery at normal room temperature, the device is ready to shock within 20 seconds if the initial rhythm finding is "Shock Advised"	
Output Energy (Biphasic): User configurable, sequence of three sequential shock levels ranging from 150-360 joules (200-360 joules, Japan)	

* Note: ±5% accuracy applies when disposable therapy electrodes are attached. Energy output is limited to the available energy which results in delivery of 360 joules into 50 ohms.	
cpfMAX technology setup options (items marked with * are default settings):	
- Stacked shocks: off, on - Initial CPR: off*, analyze first, CPR first - Preshock CPR: off*, 15, 30 seconds - Pulse check: never*, after second no shock advised, after every no shock advised, always - CPR time 1 & 2: 15, 30, 45, 60, 90, 120*, 180 seconds, 30 minutes	

PACER	
Pacing Mode: Demand pacing (user configurable) defaults (user configurable)	
Pacing Rate: 40 to 170 bpm	
Rate Accuracy: +/- 1.5%	
Output Waveform: Monophasic pulse (user configurable) current pulse (20 + 1)	
Output Current: 0 to 20 mA	
Pause: Pacing pulse for 4 when activated	
Refractory Period: 200 ms	

ENVIRONMENTAL	
Temperature, Operating SpO2: 5° to 45°C (41° to 113°F)	
Temperature, Non-operating (-4° to 140°F) except (-4° to 140°F) except	
Relative Humidity, Operating	
Atmospheric Pressure, (0 to 4572m) (0 to 15 to 4572ft)	
Water Resistance, Operating IEC 60529 (with battery)	
EMC: IEC 60601-1-2: Medical Equipment - Collateral Standard: Requirements and Tests	
IEC 60601-2-4:2002; Clause 36, Particular Cardiac Defibrillators monitors	
Shock (drop): Five drops a steel surface	
Vibration: MIL-STD-883C - category 4, and Ground Mobile -	

AC POWER ADAPTABILITY	
Function	
Dimensions: 27.7 x 16 x 10 cm (10.9 x 6.3 x 3.9 in)	
Weight: ± 2.3kg (<5 lb)	
Charge Time (with fully charged battery): 10 to 12 hours	
FASTPAK and FASTPAK NiCd: 2.1 hours	
LIFEPAK NiCd: 3.0 hours	
LIFEPAK NiCd: 3.0 hours	
LIFEPAK SLA: 6 hours	
AC Input: Accepts line voltage from 47 to 63Hz (domestic) 380 to 420Hz (military)	
Fuses: Two 250V fuses, 250mA (T2.5A) 220 to 240V: T2.5A)	

Environmental	
IPX4 per IEC 60529	
Altitude, Operating: To 10,000 ft (3,048 m)	
Altitude, Non-operating: To 10,000 ft (3,048 m)	
Humidity: 5 to 95% non-condensing	
Temperature, Operating: 5 to 150°F (-15 to 65°C)	
Temperature, Storage: (-4° to 150°F) (follows IEC 60601-1-2)	



Field-proven defibrillator/monitor in the world.



Experience the leader
that has made L1
the clear favorite
world—for every L1

For more than 50 years, Physio-Control, maker of the renowned LIFEPAK defibrillators, has been developing technologies and designing devices that are legendary among first response professionals, clinical care providers and citizens everywhere.

REFERENCES

- 1 Stiel IG, Walker RG, Nesbitt LP, et al. Biphasic Trial: A randomized comparison of fixed lower versus escalating higher energy levels for defibrillation in out-of-hospital cardiac arrest. *Circulation*. 2007;115:1511-1517.
- 2 Chapman FW, Walker RG, Koster RW. Use of 360 joule biphasic shocks for initial and recurrent ventricular fibrillation in prehospital cardiac arrest [ERC abstract O-33]. *Resuscitation*. 2006;69:49-50.
- 3 Walsh SJ, McClelland AJJ, Owen CG, et al. Efficacy of distinct energy delivery protocols comparing two biphasic defibrillators for cardiac arrest. *AM J Cardiol*. 2004;94:378-380.
- 4 McNamara RL, Wang W, Herrin J, et al. Effect of door-to-balloon time on mortality in patients with ST-segment elevation myocardial infarction. *J Am Coll Cardiol*. 2006;47:2180-2186.
- 5 Antman EM, Anbe DT, Armstrong PW, et al. ACC/AHA guidelines for the management of patients with ST-elevation myocardial infarction: executive summary: a report of the ACC/AHA Task Force on Practice Guidelines (Committee to Revise the 1999 Guidelines on the Management of Patients With Acute Myocardial Infarction). *Circulation*. 2004;110:588-636.
- 6 Bradley EH, Herrin J, Wang Y, et al. Strategies for reducing the door-to-balloon time in acute myocardial infarction. *N Engl J Med*. 2006;355:2308-2320.
- 7 Swor R, Hegerberg S, McHugh-McNally A, et al. Prehospital 12-lead ECG: efficacy or effectiveness? *Prehosp Emerg Care*. 2006;10:374-377.

To Order Call:

AEDSuperstore®

World's Largest Automated External Defibrillator Source

800.544.0048

help@aeds.com

www.AEDSuperstore.com

1800 US Hwy 51 N | Woodruff, WI 54568

**PHYSIO
CONTROL**